



PFAS – toxische Profile im Wasser



Dr. Peter Behnisch (BDS)

BioDetection Systems BV

Science Park 406

1098 XH Amsterdam

The Netherlands




Übersicht

- **Einführung zu PFAS**
- **Einführung über BIO-Analytik**
- **In vitro toxische Profile der PFAS**
- **Komplexe Mischungen von PFAS (in PFOA-Equivalenten)**
- **Erste Ergebnisse des PFAS CALUX in Wasser**
- **Zusammenfassung**


PFAS – bekannte toxische Wirkungen im Mensch

— High certainty
- - - Lower certainty







PFOS in Aal (historische Probenbank)






Mischung aus Bekannten vs. Unbekannten Stoffe im Wasser



Basis-Set: Wasser-Qualität in 5 Farben






Was können diese TOXISCHE Profile bedeuten?

Krankheit	Mechanismus	Chemikalien	CALUX Tests
Chlorakne	AhR Rezeptor	Dioxine, dl-PCBs	DR CALUX
Fettsucht	PPAR Rezeptoren	TBT, PFOA/PFOS	PPAR CALUX
Brustkrebs, Menopause	Estrogen sensitive Zellen	E2, EE2, E3, BPA, Phthalate, Pesticide	ER CALUX
Prostatakrebs	Androgen sensitive Zellen	BPA, Anabolic steroide	AR CALUX
Unfruchtbarkeit	endocrine dependent sex organ function/development	Phthalate, BPA, NP, OP,	ER, AR, TR, PR CALUX
Hypersensitivitaet	Oxidativer Stress	Pestizide	Nrf2 CALUX
Krebs Allgemein	Genotoxizitaet, DNA Schaedigung	BaP, Dioxins, PCBs	P53, PAH CALUX, Dioxin CALUX
Asthma, Allergie	Glucocorticoid sensitiv	Dexamethason, Pharmakas	GR CALUX
Thyroid gland	TTR-TR	PFAS	PFAS CALUX
Metabolisches Syndrom	complex	Dioxine, dl-PCBs, Hormone	DR und ER CALUX
Leukemie	Estrogen sensitiv	Dioxine	DR CALUX



Probenbehandlung und Bioassay Messung





Übersicht

- Einführung über BIO-Analytik
- **In vitro toxische Profile der PFAS**
- Komplexe Mischungen von PFAS in PFOA-Equivalenten
- Erste Ergebnisse des PFAS CALUX in Wasser
- Zusammenfassung


**PERFLUOROALKYL ACIDS HAVE NO ENDOCRINE DISRUPTING ACTIVITIES
VIA HUMAN STEROID HORMONE RECEPTORS AS DETERMINED BY A PANEL
OF CALUX® BIOASSAYS**

Sonneveld E¹, Man HY¹, Pieterse B¹, Brouwer A¹, Behnisch P¹, Van der Burg B¹

¹BioDetection Systems B.V. (BDS), 1098 SM, Amsterdam, the Netherlands

Compound	AR agonistic	ER α agonistic	ER β agonistic	PR agonistic	GR agonistic	AR antagonistic	ER α antagonistic	ER β antagonistic	PR antagonistic	GR antagonistic
	logEC50 (M)	logEC50 (M)	logEC50 (M)	logEC50 (M)	logEC50 (M)	logIC50 (M)	logIC50 (M)	logIC50 (M)	logIC50 (M)	logIC50 (M)
PFOS	>5	>5	>5	>5	>5	>5	>5	>5	>5	>5
PFHxA	>5	>5	>5	>5	>5	>5	>5	>5	>5	>5
PFOA	>5	>5	>5	>5	>5	>5	>5	>5	>5	>5
PFNA	>5	>5	>5	>5	>5	>5	>5	>5	>5	>5
PFHpA	>5	>5	>5	>5	>5	>5	>5	>5	>5	>5

PPAR α CALUX
(human U2OS)



Behnisch PA, Besselink H , Van der Linden S, DeVoogt P , Brouwer A.
Organohalogen Compounds 74, 186 (2012)

Chemical	Exposure levels in dust samples		Bioactivity classification (RPF)					
	Percentage of samples >MDL	Median [range (ng/g)]	Median percentage of class sum		PPARγ antagonism (Huang et al. 2016)	Thyroid hormone receptor β antagonism (Huang et al. 2016)	Androgen receptor antagonism (Huang et al. 2016)	Estrogen receptor α agonism (Huang et al. 2016)
			PPARγ antagonism (Huang et al. 2016)	Thyroid hormone receptor β antagonism (Huang et al. 2016)	Androgen receptor antagonism (Huang et al. 2016)	Estrogen receptor α agonism (Huang et al. 2016)	Thyroid hormone transport interference (Besselink 2020) ^a	
Per- and polyfluoroalkyl substances (PFAS)								
PFHxA	97.8	193 (<MDL–2,980)	66	Inactive	Inactive	Inactive	Inactive	Active (0.044) [19]
PFOS	97.8	15.2 (<MDL–296)	5.4	Active (0.55) [51]	Active (0.52) [50]	Inactive	Inactive	Active (0.85) [30]
PFOA	73.9	7.63 (<MDL–1,520)	4.5	Active (0.4) [25]	Inactive	Inactive	Inactive	Active (0.37) [11]
PFHxS	63.0	1.82 (<MDL–23.7)	<1	Unknown	Unknown	Unknown	Unknown	Active (1) [3.5]
FOSA	60.9	3.26 (<MDL–236)	1.5	Inactive	Active (0.59) [14]	Inactive	Active (0.39) [100]	Active (0.33) [3.6]
PFHpA	52.2	0.918 (0–1,760)	<1	Inactive	Inactive	Inactive	Inactive	Active (0.35) [2.0]
PFPeA	32.6	<MDL (<MDL–455)	<1	Unknown	Unknown	Unknown	Unknown	Active (0.013) [<1]
PFNA	30.4	<MDL (<MDL–1,480)	<1	Active (0.76) [2.3]	Active (0.47) [1.2]	Inactive	Inactive	Active (0.13) [<1]
PFBS	30.4	<MDL (<MDL–16.1)	<1	Unknown	Unknown	Unknown	Unknown	Active (0.028) [<1]
PFDS	10.9	<MDL (<MDL–12.5)	<1	Unknown	Unknown	Unknown	Unknown	Unknown
PFBA	4.35	<MDL (<MDL–155)	<1	Unknown	Unknown	Unknown	Unknown	Active (0.0003) [<1]
PFDA	4.35	<MDL (<MDL–35.0)	<1	Inactive	Active (0.29) [<1]	Inactive	Inactive	Active (0.033) [<1]
PFUnDA	0	<MDL (<MDL–< MDL)	<1	Active (1) [4.2]	Inactive	Inactive	Inactive	Active (0.017) [<1]
PFDoDA	0	<MDL (<MDL–< MDL)	<1	Unknown	Unknown	Unknown	Unknown	Active (0.0037) [<1]
N-MeFOSAA	0	<MDL (<MDL–< MDL)	<1	Unknown	Unknown	Unknown	Unknown	Unknown



In vitro Toxizitäts-Equivalente für PFAS (2020)

Testing of 13 PFAS standard compounds using the PFAS CALUX® bioassay.

Final report

Requested by: Dr. Tamara Grummt


Umweltbundesamt
Bismarckplatz 1
14193 Berlin
Germany

Phone: +49 3743776354
E-mail: Tamara.Grummt@uba.de




Transthyretin Rezeptor (TTR)-TR β basierter PFAS CALUX Bioassay


Step 1 Inkubation (TTR Verdringung)



Step 2 Separierung TTR-bound vs. freie T4-Stoffe



Step 3 TR β -CALUX analysis





In vitro Toxizitäts-Equivalente für PFAS (2023)

Substanz	Abkürzung	CAS Nummer	RPF (EC)	PFAS CALUX
Perfluoroctanoic acid	PFOA	335-67-1	1	1
Perfluorooctane sulfonic acid	PFOS	1763-23-1	2	2
Perfluorohexane sulfonic acid	PFHxS	355-46-4	0.6	
Perfluorononanoic acid	PFNA	375-95-1	10	
Perfluorobutanoic acid	PFBA	375-22-4	0.05	0.0018
Perfluoropentanoic acid	PFPeA	2706-90-3	0.03	0.08
Perfluorohexanoic acid	PFHxA	307-24-4	0.01	0.19
Perfluoroheptanoic acid	PFHpA	375-85-9	0.505	1.4
Perfluorodecanoic acid	PFDA	335-76-2	7	0.12
Perfluoroundecanoic acid	PFUnDA	2058-94-8	4	
Perfluorododecanoic acid	PFDoDA	307-55-1	3	
Perfluorotridecanoic acid	PFTrDA	72629-94-8	1.65	
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.3	
Perfluorohexadecanoic acid	PFHxDA	67905-19-5	0.02	
Perfluoroctadecanoic acid	PFODA	16517-11-6	0.02	
Perfluorobutane sulfonic acid	PFBS	375-73-5	0.001	0.052
Perfluoropentane sulfonic acid	PFPeS	2706-91-4	0.3005	
Perfluorohexane sulfonic acid	PFHxS		Not listed (0.6)	1.6
Perfluoroheptane sulfonic acid	PFHpS	375-92-8	1.3	1.0
Perfluorodecane sulfonic acid	PFDS	335-77-3	2	
Ammonium perfluoro (2-methyl-3-oxahexanoate)	HFPO-DA or Gen X	62037-80-3	0.06	
Propanoic Acid / Ammonium 2,2,3-trifluoro-3-(1,1,2,2,3,3-hexafluoro-3-(trifluoromethoxy)propoxy)propanoate	ADONA	958445-44-8	0.03	
6:2 Fluorotelomer sulfonic acid	6:2 FTS	27619-97-2	Not listed (0.2)	0.019
2- (Perfluorohexyl)ethyl alcohol	6:2 FTOH	647-42-7	0.02	
2-(Perfluoroctyl)ethanol	8:2 FTOH	678-39-7	0.04	
Acetic acid / 2,2-difluoro-2-((2,2,4,5-tetrafluoro-5-(trifluoromethoxy)-1,3-dioxolan-4-yl)oxy)	C6O4	1190931-41-9	0.06	

Relative Potency Factors von 24 PFAS gemäß "Proposal for a Directive of the European Parliament and of the Council amending Directive 2000/60/EC establishing a framework for Community action in the field of water policy, Directive 2006/118/EC on the protection of groundwater against pollution and deterioration and Directive 2008/105/EC on environmental quality standards in the field of water policy" in comparison to in vitro PFAS CALUX (state 01-01-2023)



PFOA Toxizitäts-Equivalente für PFAS: Trinkwasser (DE, 2017) vs Grundwasser (EU, 2023)

Bundesgesundheitsbl 2017 · 60:350–352
 DOI 10.1007/s00103-016-2508-3
 Online publiziert: 2. Januar 2017
 © Springer-Verlag Berlin Heidelberg 2017

Empfehlung des Umweltbundesamtes

Fortschreibung der vorläufigen Bewertung von per- und polyfluorinierten Chemikalien (PFC) im Trinkwasser

PFAS	TW/GOW (µg/L)	PFAS CALUX		EU RPF	Summe PFAS (EU) [µg PFOA-eq/L]
		RPF	[µg PFOA-eq/L]		
PFBA	C4	10	0.0018	0.018	0.5
PPPeA	C5	3	0.08	0.24	0.09
PFHxA	C6	6	0.19	1.14	0.06
PFHpA	C7	0.3	1.4	0.42	0.1515
PFOA	C8	0.1	1	0.1	0.1
PFNA	C9	0.06	0.32	0.0192	0.6
PFdC A	C10	0.1	0.12	0.012	0.7
PFUnDA	C11	0.1	0.052	0.0052	0.4
PFDoDA	C12	0.1	0.01	0.00098	0.3
PFTrDA	C13	0.1	0.075	0.0075	0.165
PFTeDA	C14	0.1	0.019	0.0019	0.03
PFBS	C4	6	0.052	0.31	0.006
PPPeS	C5	6		0.3005	1.803
PFHxS	C6	0.1	1.6	0.16	0.06
PFHpS	C7	0.3	1	0.3	0.39
PFOS	C8	0.1	2	0.2	0.2
PFNS	C9	0.1		2	0.2
PFDS	C10			2	0
Summe PFAS			2.9		5.8

Brussels, 26.10.2022
 COM(2022) 540 final

ANNEXES 1 to 6

ANNEXES

to the

Proposal for a Directive of the European Parliament and of the Council
 amending Directive 2000/60/EC establishing a framework for Community action in the field of water policy, Directive 2006/118/EC on the protection of groundwater against pollution and deterioration and Directive 2008/105/EC on environmental quality standards in the field of water policy

{SEC(2022) 540 final} - {SWD(2022) 540 final} - {SWD(2022) 543 final}

(1) [E ntr y] N°	(2)	(3)	(4)	(5)	(6)
	Name of substance	Category of substances	CAS number (¹)	EU number (²)	Quality Standard (³) [µg/L unless otherwise indicated]
3	Per- and poly-fluorinated alkyl substances (PFAS) - sum of 24 (⁴)	Industrial substances	See table note 6	See table note 6	0,0044 (⁷)

Faktor 1000 Differenz !



In vitro toxische Profile der PFAS (2022-2023)




PROMISCES

Preventing Recalcitrant Organic Mobile Industrial chemicals for Circular Economy in the Soil-sediment-water system

- 20 PFAS wurden auf 13 CALUX Bioassays getestet von BDS.
- Weitere 25 PFAS werden auf die 7 meist verdächtigen CALUX Bioassay getestet (Cytotox, oxidative stress Nrf2, early warning PXR, anti-TR, anti-PPAR γ , PPAR α und PFAS CALUX) von BDS
- Verschieden industriellen Produkte (GenX, ADONA, ScotchGuard), werden mit demselben CALUX Panel getestet.
- UBA testet 45 PFAS auf die biologischen Endpunkte wie z. Bsp. Zelltod, Genotoxizität (z. Bsp. Ames, Micronuclei), oxidative Stress (ROS), Hormone (ER, AR, TR) und Neurotoxizität (MEA).
- Z. Zt. werden die besten Konzentrations-Wirkungs-Kurven der meisten PFAS auf die sensitivste Manier im PFAS CALUX erfasst







Feuerlöschschäume in PFOA-eq/I

Environment International 157 (2021) 106791



Contents lists available at [ScienceDirect](#)

Environment International

journal homepage: www.elsevier.com/locate/envint



Developing potency factors for thyroid hormone disruption by PFASs using TTR-TR β CALUX® bioassay and assessment of PFASs mixtures in technical products

Peter A. Behnisch ^a, Harrie Besselink ^a, Roland Weber ^b, Wolfram Willand ^c, Jun Huang ^d, Abraham Brouwer ^{a,c,*}





Feuerlöschschäume in PFOA-eq/l



Cite This: Environ. Sci. Technol. Lett. 2019, 6, 259–264

pubs.acs.org/journal/estlcu

Letter

Per- and Polyfluoroalkyl Substances in Representative Fluorocarbon Surfactants Used in Chinese Film-Forming Foams: Levels, Profile Shift, and Environmental Implications

Mehvish Mumtaz, Yixiang Bao, Liquan Liu, Jun Huang,*^{ID} Giovanni Cagrietta,*^{ID} and Gang Yu^{ID}

State Key Joint Laboratory of Environment Simulation and Pollution Control (SKJLESPC), Beijing Key Laboratory for Emerging Organic Contaminants Control (BKLEOC), School of Environment, Tsinghua University, Beijing 100084, China

Untreated FFF																	UHPLC-MS/MS		PFAS-CALUX			
	[g/l]	4A	5A	6A	7A	8A	9A	10A	4S	5S	6S	7S	8S	9S	10S	PFOSA	42FTS	62FTS	82FTS	g PFOA-TEQ/I	Dominant	Impact 6:2 FTs
VF-368 (2013)	0.29	0.17	0.10	0.09	0.30	0.00	0.00	0.01	0.03	0.53	0.55	22.0	0.08	0.58	na	0.00	0.00	0.00	41	96 % PFOS		45
VF-570 (2013)	1.62	0.93	0.49	0.30	0.95	0.00	0.00	0.20	0.38	7.64	5.45	87.3	0.52	0.00	na	0.00	0.00	0.00	174	90 % PFOS		57
VF-9128 (2019)	0.26	0.09	1.40	0.00	0.00	0.00	0.00	0.00	0.00	0.08	0.0	0.00	0.0	0.00	na	0.00	22.4	0.00	0,70		67% 6:2 FTs	1,6
FUMETROL 21	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.0	0.00	0.0	0.00	na	0.07	1.16	0.01	0,04		61% 6:2 FTs	1,4
TEF-factors PFAS CALUX	0.0019	0.0350	0.16	1.2	1.0	0.33	0.11	0.05	na	2.0	0.0013	1.8	0.71	na	0.021	na						
na = not yet analysed																						
Heat-Top Assay treated																	UHPLC-MS/MS		PFAS-CALUX			
[g/l]	4A	5A	6A	7A	8A	9A	10A	4S	5S	6S	7S	8S	9S	10S	PFOSA	42FTS	62FTS	82FTS	g PFOA-TEQ/I	Dominant	Impact 6:2 FTs	g PFOA-BEQ/I
VF-368 (2013)	2.1	1.6	2.2	2.8	12	0.29	0.00	0.009	0.027	0.53	0.55	22	0.08	0.58	na				57	70 % PFOS		270
VF-570 (2013)	2.1	3.4	2.2	2.7	6.4	0.23	0.00	0.20	0.38	7.6	5.5	87	0.52	0.00	na				183	86 % PFOS		180
VF-9128 (2019)	1.6	2.8	0.82	0.13	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	na				0,4	42% PFHpA		<LOQ
FUMETROL 21	4.9	12	3.4	0.63	0.02	0.00	0.00	0.00	0.01	0.00	0.01	0.00	0.00	0.00	na	0.00	0.00	0.00	1,8	39% PFHpA	0% 6:2 FTs	<LOQ

Behnisch et al 2021. Environmental International



Oberflächen- und Abfall-Wasser in PFOA-eq/l (SETAC, IWA und DIOXIN 2022)

Evaluation of thyroid hormone disruption by PFAS in WWTP influent/effluent and surface waters in the Netherlands

Peter A. Behnisch¹, Harrie Besselink¹, Dorien ten Hulscher², Anne Jans², Carmen Hogendoorn², John Hin², Carolien van der Wielen¹, Abraham Brouwer¹


¹BioDetection Systems BV, Science Park 406, 1098 XH, Amsterdam, the Netherlands,

²Rijkswaterstaat Water, Verkeer en Leefomgeving (RWS-WVL), Griffioenlaan 2, 3526 LA, Utrecht, the Netherlands.

H. Besselink¹, P. Behnisch¹, D. ten Hulscher², A. Jans², C. Hogendoorn², J. Hin², C. van der Wielen¹, A. Brouwer¹


¹ BioDetection Systems BV, Science Park 406, 1098 XH, Amsterdam, the Netherlands,

² Rijkswaterstaat Water, Verkeer en Leefomgeving (RWS-WVL), Griffioenlaan 2, 3526 LA, Utrecht, the Netherlands.





Oberflächenwasser (NL): PFAS CALUX (Sum in vitro PFAS) vs LC/MS (13 PFAS)



	PFAS CALUX (ug total PFOA eq./l)	LC/MS (PFAS 13) (ug PFOA eq./l)
5% Percentile	0.56	0.0038
80% Percentile	4.3	0.035



Scientific Committee on Health, Environmental and Emerging Risks
SCHEER

Scientific Opinion on "Draft Environmental Quality Standards for Priority Substances under the Water Framework Directive"

PFAS




The SCHEER adopted this document
via written procedure on 18 August 2022.



Oberflächen und Abfall-Wasser
PFAS CALUX (Summe in vitro PFAS) vs LCMS (13 PFAS)


Rot = Oberflächenwasser
Blau = Kläranlagenabwasser





Oberflächen und Abfall-Wasser PFAS CALUX (Sum in vitro PFAS) vs LCMS (13 PFAS)

Summe PFOA äquivalente
Lila = LC/MS/MS
Grün = PFAS CALUX



Percentile (%)	PFAS CALUX (ug PFOA eq./l water)	LCMS (PFAS 13; upperbound) (ug PFOA eq./l water)
5	0.71	0.0094
80	22	0.19

**PFOA-Äquivalente für
Oberflächenwasser und
Abwasser:**

Um diese PFAS-EQ-Gesamtlücken zu füllen, müssen mehr PFAS durch chemische Analyse analysiert und mehr In-vitro-Potenzfaktoren für diese PFAS-Verbindungen ermittelt werden.



In der Kürze liegt die Würze ...warum toxische Profile in Wasser ?

- **Wirkungs- und Summen-bezogen**
- **Ethisch (kein Tierversuch, NAM)**
- **Relevant für Mensch (Menschliche Zellen)**
- **Luminometer-basiert**
- **Relevant für Green Deal & Zero-Pollution & Toxic-freien Strategien**





Herzlichen Dank für Ihr freundliches Interesse!



BDS
BioDetection Systems BV

Science Park 406
1098 XH Amsterdam
The Netherlands

Tel.: +31 (0)20 4350750
E-mail: Behnisch@bds.nl